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## Theoretical Background

1. History .....	1
1.1 Known for More Than 100 Years: the Reflux Condenser .....	1
1.2 Reflux Heating – the Perfect Technology? .....	1
2. How Can We Accelerate Chemical Transformations? .....	1
2.1 Activation Energy and the Arrhenius Law .....	1
3. Superheating in Closed Vessels .....	3
3.1 Intelligent Process Control and Safety .....	3

## The Monowave 50 Synthesis Reactor

1. What is Monowave 50? .....	4
2. Specifications .....	5
3. Handling: It's That Simple! .....	5

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## Anhang

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## 1. History

### 1.1 Known for More Than 100 Years: the Reflux Condenser

Among the first things that chemistry students notice about organic chemistry is that many organic reactions would proceed rather slowly without additional heating. What further complicates things is that heating of volatile, highly flammable organic solvents in open vessels does not only generate a safety risk, but will also very quickly lead to dry-off of the reaction batch. The historical solution to this problem still in use today is the reflux heating methodology utilizing water-cooled glass condensers to prevent solvent escaping into the open atmosphere.

### 1.2 Reflux Heating – the Perfect Technology?

Reflux condensers belong to the standard inventory of most labs and can still be considered the synthetic “workhorse”.

However, the time in educational institutions is usually limited, which makes it difficult - sometimes even impossible – for students to perform common chemical reactions that usually take several hours or even days. Furthermore, environmental issues have to be considered, since reflux experiments waste a lot of water and also chemicals due to commonly used comparably large reaction flasks.


A particular weakness of the widely used reflux heating methodology is undoubtedly the accessible temperature range (being dependent on the solvent's boiling point). Thereby, the solvent choice also co-determines the reaction temperature. A mostly unavoidable consequence of using the reflux heating method is that very long reaction times of between hours and days have to be accepted.

Furthermore, reflux experiments bear some additional risks. Pressure fluctuations in the cooling water supply and/or improperly fixed equipment (an issue encountered especially when students are involved) could lead to water leakage.

## 2. How Can We Accelerate Chemical Transformations?

### 2.1 Activation Energy and the Arrhenius Law

Figure 1.1 shows an energy diagram of a chemical reaction. **G** is the free energy of a given substance. If the free energy of the products is lower than the free energy of the educts, the reaction is energetically favored meaning that this reaction will take place spontaneously if the activation energy barrier can be overcome. **E<sub>a</sub>**, the activation energy, can be seen as an energetic barrier which needs to be overcome in order for the reaction to proceed. Therefore, some chemical reactions need to be heated in order to take place.

## Theoretical Background

There is a relationship between the reaction rate constant **k** and the reaction temperature. The Arrhenius Law describes this relationship.

$$k = A \times e^{-E_a/RT}$$

**k** = rate constant (quantifies the speed of a chemical reaction)

**A** = pre-exponential factor

**E<sub>a</sub>** = activation energy

**R** = universal gas constant

**T** = reaction temperature

**E<sub>a</sub>** and **k** are different for every chemical reaction and this is the reason why some reactions already take place at room temperature or even temperatures below 0 °C whereas others need to be heated to high temperatures. Reaction times can be calculated as a function of the reaction temperature using a simplified version of the Arrhenius equation.

### Keep in mind:

Rule of Thumb: 10 °C temperature increase = 2-fold rate acceleration

### Example:

If a reaction would take 8 hours at 80 °C, heating to 90 °C would reduce the reaction time by 50 % to 4 hours. Upon further heating, the reaction time can be further reduced, so that at 140 °C, it takes only 8 minutes to complete this reaction.

Temperature	80 °C	90 °C	100 °C	110 °C	120 °C	130 °C	140 °C
Time	8 h	4 h	2 h	1 h	30 min	15 min	8 min

### 3. Superheating in Closed Vessels

Why and how are we able to shorten the reaction times and increase our efficiency using Monowave 50? In fact, Monowave 50 works under sealed vessel conditions. If the reaction vessels are sealed, the reaction mixture can reach temperatures far above the boiling point of the used solvent(s), which leads to dramatically reduced reaction times (from several hours down to a few minutes!).

According to the Arrhenius law the simplest way of increasing the rate of a reaction is to increase the temperature, a strategy that is usually limited by the boiling point of the used solvent. If, however, a reaction mixture is heated inside a sealed reaction vessel, no solvent vapor can escape and theoretically, any reaction temperature can be realized. This method is often referred to as "superheating" of a solvent and requires the use of reinforced reaction vessels because of the built-up pressure. In a closed vessel, the pressure of the trapped liquid heated to a temperature above the boiling point equals the vapor pressure.

#### Keep in mind:

Only closed vessel systems provide reaction enhancement due to superheating above the boiling point of the used solvent(s).

If a solvent is heated above its boiling point, pressures higher than atmospheric pressure will be reached and this requires closed vessels.

As an example, water is heated to about 120 °C inside a domestic pressure cooker used in many households to achieve a faster cooking process. At this temperature, the vapor pressure of water of around 2 bar is acting on the vessel walls. Because of this pressure buildup, reactors which allow high pressures need to be specifically equipped to fulfill the necessary safety requirements.

#### 3.1 Intelligent Process Control and Safety



Monowave 50 is able to perform all necessary monitoring and control steps automatically. The reactions are not only temperature-controlled and stopped automatically after a set time, they are also automatically terminated in case of e.g. sudden pressure drop or rise (vessel over-pressurization, burst or leak) and other incidents which can be handled by the instrument without the user needing to be present. While the experiment is running, the safety lock of Monowave 50 prevents from opening the instrument cover. Whereas a stirred reflux set-up immersed in an oil bath would continue to heat and deliver water to a broken reflux setup, Monowave 50 would automatically stop heating, immediately cool the reactor cavity and display a description of the incident on the screen.

# The Monowave 50 Synthesis Reactor

## 1. What is Monowave 50?

Monowave™ 50 is the first conventional heating synthesis reactor made by Anton Paar and is especially designed for educational labs and standard experimental chemistry.

The reactor is able to heat any kind of solvents and to perform chemical reactions up to 250 °C and 20 bar. A contact sensor placed at the bottom of the cavity allows accurate temperature control during the whole experiment.

The 10 mL borosilicate glass vials are easily closed with silicone caps. After placing the vial into the reactor cavity, the system is hermetically sealed and the experiment can be started in a few steps via its intuitive capacitive touchscreen.

When a run is finished and the vial has been cooled down to a safe temperature, the cover is opened manually: while opening, any pressure generated during the reaction is safely released into the exhaust channel of Monowave. Hence, the handling risk for the user is minimized. After the system has been opened the processed reaction mixture is easily and safely accessible and no handling of pressurized vials is required.

The integrated overpressure controller allows reactions up to 20 bar to be processed safely. That means that if the pressure limit is reached during the reaction the power is reduced to maintain the reaction below the pressure limit.

Stirring efficiency, a crucial feature in chemical synthesis, is ensured by the installed magnetic stirrer which allows stirring speeds between 300 rpm and 1200 rpm.

Simplicity has been considered for Monowave 50 even in terms of cleaning and maintenance. Only vials and caps have to be cleaned after the reaction. This simplifies the laboratory work and saves time. In the unlikely event of a vessel rupture, cleaning can be performed straightforwardly.

# The Monowave 50 Synthesis Reactor

## 2. Specifications

Max. filling volume	6 mL
Max. operation pressure	20 bar
Max. temperature	250°C
Max. power	315 W
Vial material	Borosilicate glass
Cap material	Silicone
Seal material	PTFE

## 3. Handling: It's That Simple!

### Keep in mind:

Before using Monowave 50, read the Monowave 50 Instruction Manual carefully and observe all mentioned safety instructions!

Typically, doing an experiment in Monowave 50 involves the following steps:

1. Fill the G10 vial and close it with a silicone septum
2. Introduce the vial into the reactor cavity
3. Close the reactor
4. Select or create a method or program the experiment in the Quick Settings

A Monowave 50 experiment consists of **three phases**: the **heating phase** (as fast as possible (AFAP)) or with a defined time to heat to the target temperature), the **hold phase** (the target temperature is maintained for a given time) and the **cool down phase** (to reach a safe handling temperature). The latter one does automatically take place as soon as the set heating time of an experiment has elapsed.

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